

For Research Use Only

SC4MOL Polyclonal antibody, PBS Only

Catalog Number: 31159-1-PBS



Basic Information

Catalog Number:

31159-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34758

GenBank Accession Number:

BC107879

GeneID (NCBI):

6307

UNIPROT ID:

Q15800

Full Name:

sterol-C4-methyl oxidase-like

Observed MW:

34-36 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SC4MOL(sterol C4-methyl oxidase-like), an intermediate enzyme in the cholesterol biosynthetic pathway, was included in the library based on its representation in a high confidence GEO transcriptional profiling dataset as a transcript that rapidly undergoes significant expression change in response to stimulation or inhibition of EGFR. Among the validated hits that increased cell killing by common EGFR antagonists, SC4MOL was consistently one of the most effective modulators. (PMID: 23125191)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31159-1-AP (SC4MOL antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31159-1-PBS in a different storage buffer formulation.