

For Research Use Only

GANC Polyclonal antibody, PBS Only

Catalog Number: 31163-1-PBS



Basic Information

Catalog Number:

31163-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35067

GenBank Accession Number:

BC059406

GeneID (NCBI):

2595

UNIPROT ID:

Q8TET4

Full Name:

glucosidase, alpha; neutral C

Calculated MW:

104 kDa

Observed MW:

100-105 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

GANC (Glucosidase Alpha, Neutral C) is a Protein Coding gene. Diseases associated with GANC include Pompe Disease and Autosomal Recessive Limb-Girdle Muscular Dystrophy Type 2A.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31163-1-AP (GANC antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31163-1-PBS in a different storage buffer formulation.