

For Research Use Only

UBAP1L Polyclonal antibody, PBS Only

Catalog Number: 31189-1-PBS



Basic Information

Catalog Number:

31189-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34816

GenBank Accession Number:

NM_001163692

GeneID (NCBI):

390595

UNIPROT ID:

F5GYI3

Full Name:

similar to ubiquitin-associated protein 1 (predicted)

Calculated MW:

41 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Caco-2 cells were subjected to SDS PAGE followed by western blot with 31189-1-AP (UBAP1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31189-1-PBS in a different storage buffer formulation.