

For Research Use Only

TRIO Polyclonal antibody

Catalog Number: 31192-1-AP



Basic Information

Catalog Number:

31192-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34954

GenBank Accession Number:

NM_007118

GeneID (NCBI):

7204

UNIPROT ID:

O75962

Full Name:

triple functional domain (PTPRF interacting)

Calculated MW:

347 kDa

Observed MW:

347 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31192-1-AP (TRIO antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.