

For Research Use Only

CLK4 Polyclonal antibody

Catalog Number: 31205-1-AP



Basic Information

Catalog Number:

31205-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34446

GenBank Accession Number:

NM_020666

GeneID (NCBI):

57396

UNIPROT ID:

Q9HAZ1

Full Name:

CDC-like kinase 4

Calculated MW:

57 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : EC109 cells, K-562 cells, mouse testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31205-1-AP (CLK4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.