

For Research Use Only

SLC24A2 Polyclonal antibody

Catalog Number: 31289-1-AP



Basic Information

Catalog Number:

31289-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35294

GenBank Accession Number:

NM_020344

GeneID (NCBI):

25769

UNIPROT ID:

Q9UI40

Full Name:

solute carrier family 24 (sodium/potassium/calcium exchanger), member 2

Calculated MW:

661 aa, 74 kDa

Observed MW:

73-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : ARPE-19 cells, hTERT-RPE1 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31289-1-AP (SLC24A2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.