

For Research Use Only

WDR21A Polyclonal antibody

Catalog Number: **31440-1-AP**



Basic Information

Catalog Number:

31440-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35426

GenBank Accession Number:

BC018979

GeneID (NCBI):

26094

UNIPROT ID:

Q8WV16

Full Name:

WD repeat domain 21A

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HeLa cells, Jurkat cells, U-251 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31440-1-AP (WDR21A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.