

For Research Use Only

HERC1 Polyclonal antibody, PBS Only

Catalog Number: 31535-1-PBS



Basic Information

Catalog Number:

31535-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33624

GenBank Accession Number:

NM_003922

GeneID (NCBI):

8925

UNIPROT ID:

Q15751

Full Name:

hect (homologous to the E6-AP
(UBE3A) carboxyl terminus) domain
and RCC1 (CHC1)-like domain (RLD) 1

Calculated MW:

532 kDa

Observed MW:

532-550 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 31535-1-AP (HERC1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31535-1-PBS in a different storage buffer formulation.