

For Research Use Only

MXRA8 Polyclonal antibody

Catalog Number: 31581-1-AP



Basic Information

Catalog Number: 31581-1-AP	GenBank Accession Number: BC017312	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 300 ug/ml by Nanodrop;	GeneID (NCBI): 54587	Recommended Dilutions: WB: 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9BRK3	
Isotype: IgG	Full Name: matrix-remodelling associated 8	
Immunogen Catalog Number: AG33429	Calculated MW: 49 kDa	
	Observed MW: 70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COS-7 cells, mouse ovary tissue, rat ovary tissue
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31581-1-AP (MXRA8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.