

For Research Use Only

JPH2 Polyclonal antibody

Catalog Number: 31589-1-AP



Basic Information

Catalog Number:

31589-1-AP

Size:

150ul , Concentration: 480 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36314

GenBank Accession Number:

BC043206

GeneID (NCBI):

57158

UNIPROT ID:

Q9BR39

Full Name:

junctophilin 2

Observed MW:

74-97 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Caco-2 cells, mouse skeletal muscle tissue, rat skeletal muscle tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31589-1-AP (JPH2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.