

For Research Use Only

# RGPD5 Polyclonal antibody, PBS Only

Catalog Number: 31636-1-PBS



## Basic Information

**Catalog Number:**

31636-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG36253

**GenBank Accession Number:**

NM\_001037866

**GeneID (NCBI):**

84220

**UNIPROT ID:**

Q99666

**Full Name:**

RANBP2-like and GRIP domain containing 5

**Calculated MW:**

198kDa

**Purification Method:**

Antigen affinity Purification

## Applications

**Tested Applications:**

IHC, Indirect ELISA

**Species Specificity:**

human, mouse

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 31636-1-AP (RGPD5 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31636-1-PBS in a different storage buffer formulation.