

For Research Use Only

NUAK2 Polyclonal antibody

Catalog Number: 31718-1-AP



Basic Information

Catalog Number: 31718-1-AP	GenBank Accession Number: NM_030952	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 380 ug/ml by Nanodrop;	GeneID (NCBI): 81788	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q9H093	
Isotype: IgG	Full Name: NUAK family, SNF1-like kinase, 2	
Immunogen Catalog Number: AG35302	Calculated MW: 70kDa	
	Observed MW: 70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HEK-293T cells, HeLa cells
Species Specificity: human	

Background Information

NUAK2, also known as OMPK2 and SNARK, is the fourth member of the AMPK kinase family. The molecular mass of NUAK2 is 69 kDa and it plays a key role in neural tube closure during embryonic development (PMID: 32845958).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31718-1-AP (NUAK2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.