

For Research Use Only

ZFP42 Polyclonal antibody

Catalog Number: 31790-1-AP



Basic Information

Catalog Number: 31790-1-AP	GenBank Accession Number: BC160056	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 500 ug/ml by Nanodrop;	GeneID (NCBI): 132625	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q96MM3	
Isotype: IgG	Full Name: zinc finger protein 42 homolog (mouse)	
Immunogen Catalog Number: AG36112	Observed MW: 47 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HL-60 cells, LNCaP cells, PC-3 cells, Raji cells
Species Specificity: human	

Background Information

The ZFP42, also known as REX1, is a transcriptional regulator and a member of the zinc finger protein family. It is specifically highly expressed in pluripotent stem cells, such as embryonic stem cells, and is one of the key molecules for maintaining cell pluripotency and self-renewal capacity. Its expression level rapidly declines with cell differentiation, making it widely used as a reliable marker for the pluripotent state. Studies have shown that ZFP42 is involved in inhibiting differentiation programs and maintaining stemness by regulating the transcription of downstream target genes. Additionally, its abnormal expression in various tumors suggests that it may be involved in tumorigenesis and development, making it an important subject in the research of regenerative medicine, developmental biology, and oncology.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31790-1-AP (ZFP42 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.