

For Research Use Only

PSMD1 Polyclonal antibody

Catalog Number: 31912-1-AP



Basic Information

Catalog Number:

31912-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36142

GenBank Accession Number:

BC094720

GeneID (NCBI):

5707

UNIPROT ID:

Q99460

Full Name:

proteasome (prosome, macropain)
26S subunit, non-ATPase, 1

Calculated MW:

953 aa, 106 kDa

Observed MW:

106 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : A431 cells, Jurkat cells, SKOV-3 cells, mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31912-1-AP (PSMD1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.