

For Research Use Only

PITHD1 Polyclonal antibody

Catalog Number:31973-1-AP



Basic Information

Catalog Number:

31973-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37120

GenBank Accession Number:

BC017208

GeneID (NCBI):

57095

UNIPROT ID:

Q9GZP4

Full Name:

chromosome 1 open reading frame

128

Calculated MW:

24 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HL-60 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31973-1-AP (PITHD1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.