

For Research Use Only

CD44 Polyclonal antibody

Catalog Number: 31981-1-AP



Basic Information

| | | |
|---|--|---|
| Catalog Number: 31981-1-AP | GenBank Accession Number: XM_006234631.3 | Purification Method: Antigen affinity Purification |
| Size: 150ul , Concentration: 500 ug/ml by Nanodrop; | GeneID (NCBI): 25406 | Recommended Dilutions: WB: 1:500-1:2000 IF-P: 1:50-1:500 |
| Source: Rabbit | UNIPROT ID: O70509 | |
| Isotype: IgG | Full Name: Cd44 molecule | |
| Immunogen Catalog Number: EG1458 | Calculated MW: 40 kDa | |
| | Observed MW: 80-95 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB, IF-P, ELISA | Positive Controls: |
| Species Specificity: mouse, rat | WB : mouse spleen tissue, rat lung tissue, rat spleen tissue IF-P : rat spleen tissue, |

Background Information

CD44 is a type I transmembrane glycoprotein expressed on embryonic stem cells and in various levels on other cell types including connective tissues and bone marrow. CD44 expression is also upregulated in subpopulations of cancer cells and is recognized as a molecular marker for cancer stem cells (PMID: 29747682). It is a cell-surface receptor that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands (PMID: 10694938). Adhesion with HA plays an important role in cell migration, tumor growth and progression. CD44 is also involved in lymphocyte activation, recirculation and homing, and in hematopoiesis.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31981-1-AP (CD44 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat spleen tissue using CD44 antibody (31981-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat spleen tissue using CD44 antibody (31981-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).