

For Research Use Only

LANCL3 Polyclonal antibody, PBS Only

Catalog Number:32087-1-PBS



Basic Information

Catalog Number:

32087-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34968

GenBank Accession Number:

BC093667

GeneID (NCBI):

347404

UNIPROT ID:

Q6ZV70

Full Name:

LanC lantibiotic synthetase component C-like 3 (bacterial)

Calculated MW:

46 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32087-1-AP (LANCL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32087-1-PBS in a different storage buffer formulation.