

For Research Use Only

PRHOXNB Polyclonal antibody, PBS Only

Catalog Number:32118-1-PBS



Basic Information

Catalog Number: 32118-1-PBS	GenBank Accession Number: NM_001105577	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 646625	
Source: Rabbit	UNIPROT ID: A6NGE7	
Isotype: IgG	Full Name: parahox cluster neighbor	
Immunogen Catalog Number: AG37232	Calculated MW: 19kDa	
	Observed MW: 20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

PRHOXNB may Catalyzes the stereoselective decarboxylation of 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline (OHCU) to (S)-allantoin. PRHOXNB is expressed at extremely low levels.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32118-1-AP (PRHOXNB antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32118-1-PBS in a different storage buffer formulation.