

For Research Use Only

# OSBPL6 Polyclonal antibody

Catalog Number:32158-1-AP



## Basic Information

<b>Catalog Number:</b> 32158-1-AP	<b>GenBank Accession Number:</b> BC052259	<b>Purification Method:</b> Antigen affinity Purification
<b>Size:</b> 150ul , Concentration: 570 ug/ml by Nanodrop;	<b>GeneID (NCBI):</b> 114880	<b>Recommended Dilutions:</b> WB: 1:500-1:3000 IHC: 1:50-1:500
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9BZF3	
<b>Isotype:</b> IgG	<b>Full Name:</b> oxysterol binding protein-like 6	
<b>Immunogen Catalog Number:</b> AG37534	<b>Calculated MW:</b> 934 aa, 106 kDa	
	<b>Observed MW:</b> 106 kDa	

## Applications

### Tested Applications:

WB, IHC, ELISA

### Species Specificity:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

### Positive Controls:

WB : HEL cells,

IHC : mouse skeletal muscle tissue, mouse brain tissue, mouse lung tissue

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 32158-1-AP (OSBPL6 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 32158-1-AP (OSBPL6 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 32158-1-AP (OSBPL6 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 32158-1-AP (OSBPL6 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).