

For Research Use Only

FAM161A Polyclonal antibody

Catalog Number: 32256-1-AP



Basic Information

Catalog Number:

32256-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35562

GenBank Accession Number:

BC107163

GeneID (NCBI):

84140

UNIPROT ID:

Q3B820

Full Name:

family with sequence similarity 161, member A

Calculated MW:

716 aa, 83 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : human placenta tissue, rat liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32256-1-AP (FAM161A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.