

For Research Use Only

FAM161A Polyclonal antibody, PBS Only

Catalog Number:32256-1-PBS



Basic Information

Catalog Number:

32256-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35562

GenBank Accession Number:

BC107163

GeneID (NCBI):

84140

UNIPROT ID:

Q3B820

Full Name:

family with sequence similarity 161, member A

Calculated MW:

716 aa, 83 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32256-1-AP (FAM161A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32256-1-PBS in a different storage buffer formulation.