

For Research Use Only

SLC35F4 Polyclonal antibody, PBS Only

Catalog Number: 32258-1-PBS

Featured Product



Basic Information

Catalog Number:

32258-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34518

GenBank Accession Number:

NM_001206920

GeneID (NCBI):

341880

UNIPROT ID:

A4IF30

Full Name:

solute carrier family 35, member F4

Calculated MW:

521 aa, 58 kDa

Observed MW:

42 kDa, 65 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32258-1-AP (SLC35F4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32258-1-PBS in a different storage buffer formulation.



WB result of SLC35F4 antibody (32258-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-SLC35F4 transfected A431 cells. This data was developed using the same antibody clone with 32258-1-PBS in a different storage buffer formulation.