

For Research Use Only

LINGO3 Polyclonal antibody

Catalog Number:32267-1-AP



Basic Information

Catalog Number:

32267-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37289

GenBank Accession Number:

NM_001101391

GeneID (NCBI):

645191

UNIPROT ID:

POC6S8

Full Name:

leucine rich repeat and Ig domain containing 3

Calculated MW:

64kDa

Observed MW:

67 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A549 cells, Calu-1 cells, U-251 cells, U-87 MG cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32267-1-AP (LINGO3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.