

For Research Use Only

ZNF530 Polyclonal antibody, PBS Only

Catalog Number: 32288-1-PBS



Basic Information

Catalog Number:

32288-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37782

GenBank Accession Number:

BC060865

GeneID (NCBI):

348327

UNIPROT ID:

Q6P9A1

Full Name:

zinc finger protein 530

Calculated MW:

599 aa, 69 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

ZNF530, also known as zinc finger protein 530, is a member of the zinc finger protein family. ZNF530 is expressed in various tissues, including the brain, heart, liver, and pancreas. It has been shown to play a role in the development and progression of cancer, as well as in the regulation of normal cell development and function.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32288-1-AP (ZNF530 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32288-1-PBS in a different storage buffer formulation.