

For Research Use Only

HEATR5B Polyclonal antibody

Catalog Number: 32316-1-AP



Basic Information

Catalog Number:

32316-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36352

GenBank Accession Number:

BC150605

GeneID (NCBI):

54497

UNIPROT ID:

Q9P2D3

Full Name:

HEAT repeat containing 5B

Observed MW:

~210 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A549 cells, U-251 cells, RT-4 cells

Background Information

HEATR5B (HEAT Repeat Containing 5B) is a protein containing HEAT repeat that binds directly to the dynein tail and dynactin and has an evolutionarily conserved function in promoting the motility of adaptor protein-1 (AP1)-associated endosomal membranes (PMID: 37872872).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32316-1-AP (HEATR5B antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.