

For Research Use Only

# USP9Y Polyclonal antibody

Catalog Number: 32373-1-AP



## Basic Information

Catalog Number:

32373-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35809

GenBank Accession Number:

NM\_004654

GeneID (NCBI):

8287

UNIPROT ID:

O00507

Full Name:

ubiquitin specific peptidase 9, Y-linked

Calculated MW:

291kDa

Observed MW:

250 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : mouse brain tissue,

## Background Information

USP9Y (Ubiquitin-Specific Peptidase 9 on the Y chromosome) is a gene located in the male-specific region of the Y chromosome and is broadly expressed in male adult tissues. It encodes a functional deubiquitinating enzyme that plays a role in protein turnover (PMID: 39742627).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32373-1-AP (USP9Y antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.