

For Research Use Only

USP9Y Polyclonal antibody, PBS Only

Catalog Number: 32373-1-PBS



Basic Information

Catalog Number:

32373-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35809

GenBank Accession Number:

NM_004654

GeneID (NCBI):

8287

UNIPROT ID:

O00507

Full Name:

ubiquitin specific peptidase 9, Y-linked

Calculated MW:

291kDa

Observed MW:

250 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

USP9Y (Ubiquitin-Specific Peptidase 9 on the Y chromosome) is a gene located in the male-specific region of the Y chromosome and is broadly expressed in male adult tissues. It encodes a functional deubiquitinating enzyme that plays a role in protein turnover (PMID: 39742627).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32373-1-AP (USP9Y antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32373-1-PBS in a different storage buffer formulation.