

For Research Use Only

# ELA3B Polyclonal antibody, PBS Only

Catalog Number: 32376-1-PBS



## Basic Information

Catalog Number:

32376-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36965

GenBank Accession Number:

BC005216

GeneID (NCBI):

23436

UNIPROT ID:

P08861

Full Name:

elastase 3B, pancreatic

Calculated MW:

29 kDa

Observed MW:

29-32 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Human proelastases comprise proelastase 2A (ELA2A), proelastase 3A (ELA3A) and proelastase 3B (ELA3B).

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32376-1-AP (ELA3B antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32376-1-PBS in a different storage buffer formulation.