

For Research Use Only

CDYL2 Polyclonal antibody

Catalog Number: 32433-1-AP



Basic Information

Catalog Number: 32433-1-AP	GenBank Accession Number: BC069440	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 400 ug/ml by Nanodrop;	GeneID (NCBI): 124359	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8N8U2	
Isotype: IgG	Full Name: chromodomain protein, Y-like 2	
Immunogen Catalog Number: AG37694	Calculated MW: 506 aa, 57 kDa	
	Observed MW: 55-57 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, MCF-7 cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32433-1-AP (CDYL2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.