

For Research Use Only

ZCCHC4 Polyclonal antibody, PBS Only

Catalog Number: 32481-1-PBS



Basic Information

Catalog Number: 32481-1-PBS	GenBank Accession Number: BC157834	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 29063	
Source: Rabbit	UNIPROT ID: Q9H5U6	
Isotype: IgG	Full Name: zinc finger, CCHC domain containing 4	
Immunogen Catalog Number: AG34661	Observed MW: 55 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32481-1-AP (ZCCHC4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32481-1-PBS in a different storage buffer formulation.