

For Research Use Only

SVIP Polyclonal antibody, PBS Only

Catalog Number: 32511-1-PBS



Basic Information

Catalog Number:

32511-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36618

GenBank Accession Number:

NM_001320340

GeneID (NCBI):

258010

UNIPROT ID:

Q8NHG7

Full Name:

small VCP/p97-interacting protein

Calculated MW:

8kDa

Observed MW:

12 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

SVIP (Small VCP-interacting protein) is a 9-Kda adaptor endoplasmic reticulum protein that could bind directly to p97/VCP. SVIP is located on the ER membrane's cytosolic surface. It inhibits the ERAD pathway, which is known independently of ubiquitin, by interacting with p97/VCP, which has a central role in this pathway. Especially for the ERAD pathway, the SVIP protein has been identified as a negative regulator. SVIP was identified as a novel adapter protein for p97/VCP by observing its overexpression in cells, which causes extensive vacuolation and deformation of the ER and microtubules (PMID: 39473740).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32511-1-AP (SVIP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32511-1-PBS in a different storage buffer formulation.