

For Research Use Only

FAM40A Polyclonal antibody, PBS Only

Catalog Number:32554-1-PBS



Basic Information

Catalog Number:

32554-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38142

GenBank Accession Number:

BC121793

GeneID (NCBI):

85369

UNIPROT ID:

Q5VSL9

Full Name:

family with sequence similarity 40, member A

Calculated MW:

837 aa, 96 kDa

Observed MW:

95-110 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

FAM40A plays a role in the regulation of cell morphology and cytoskeletal organization. It is required in the cortical actin filament dynamics and cell shape.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32554-1-AP (FAM40A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32554-1-PBS in a different storage buffer formulation.