

For Research Use Only

LRIT2 Polyclonal antibody

Catalog Number: 32565-1-AP



Basic Information

Catalog Number:

32565-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38495

GenBank Accession Number:

BC144475

GeneID (NCBI):

340745

UNIPROT ID:

A6NDA9

Full Name:

leucine-rich repeat, immunoglobulin-like and transmembrane domains 2

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse retina tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32565-1-AP (LRIT2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.