

For Research Use Only

ZNF444 Polyclonal antibody, PBS Only

Catalog Number: 32588-1-PBS



Basic Information

Catalog Number:

32588-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38517

GenBank Accession Number:

NM_018337.3

GeneID (NCBI):

55311

UNIPROT ID:

Q8NOY2

Full Name:

zinc finger protein 444

Calculated MW:

35kDa 327aa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

ZNF444 (Zinc Finger Protein 444) is a member of the zinc finger protein family, which are characterized by their ability to bind DNA, RNA, and proteins through zinc finger motifs. ZNF444 is encoded by the ZNF444 gene in humans and is involved in various cellular processes, including transcriptional regulation and immune responses.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32588-1-AP (ZNF444 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32588-1-PBS in a different storage buffer formulation.