

For Research Use Only

TAF7L Polyclonal antibody, PBS Only

Catalog Number: 32590-1-PBS



Basic Information

Catalog Number:

32590-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38542

GenBank Accession Number:

NM_024885.3

GeneID (NCBI):

54457

UNIPROT ID:

Q5H9L4

Full Name:

TAF7-like RNA polymerase II, TATA box binding protein (TBP)-associated factor, 50kDa

Calculated MW:

53kDa 462aa

Observed MW:

53 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32590-1-AP (TAF7L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32590-1-PBS in a different storage buffer formulation.