

For Research Use Only

SEC14L2 Polyclonal antibody, PBS Only

Catalog Number: 32593-1-PBS



Basic Information

Catalog Number: 32593-1-PBS	GenBank Accession Number: BC058915	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 23541	
Source: Rabbit	UNIPROT ID: O76054	
Isotype: IgG	Full Name: SEC14-like 2 (S. cerevisiae)	
Immunogen Catalog Number: AG35592	Calculated MW: 46 kDa	
	Observed MW: 44 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

SEC14L2 (SEC14-like lipid-binding protein 2) is a cytoplasmic protein that belongs to the lipid-binding protein family. It plays an important role in a variety of cellular processes, including endosome division and lipid metabolism. SEC14L2 is able to bind phosphatidylinositol (PtdIns) and regulate its intracellular distribution, which is essential for normal endosome division. It has been shown that SEC14L2 regulates the endosome division process by promoting the conversion of PtdIns4P to PtdIns3P. In addition, SEC14L2 plays a role in the cholesterol biosynthesis pathway by stimulating squalene monooxygenase activity. In terms of viral infection, SEC14L2 is essential for hepatitis C virus (HCV) replication, and its expression promotes HCV replication in cells.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32593-1-AP (SEC14L2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32593-1-PBS in a different storage buffer formulation.