

For Research Use Only

HSF2BP Polyclonal antibody, PBS Only

Catalog Number:32594-1-PBS



Basic Information

Catalog Number:

32594-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37637

GenBank Accession Number:

BC000153

GeneID (NCBI):

11077

UNIPROT ID:

O75031

Full Name:

heat shock transcription factor 2 binding protein

Calculated MW:

38 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Heat shock factor 2-binding protein (HSF2BP, also known as MEILB2) is a binding protein of HSF2, which is involved in modulating HSF2 activation in testis (PMID: 9651507). Naturally occurring elevated production of HSF2BP in tumors may be a source of cancer-promoting genomic instability and a vulnerable target (PMID: 31960047; 34790765).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32594-1-AP (HSF2BP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32594-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using HSF2BP antibody (32594-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 32594-1-PBS in a different storage buffer formulation.