

For Research Use Only

HSD3B1 Polyclonal antibody

Catalog Number: 32642-1-AP



Basic Information

Catalog Number: 32642-1-AP	GenBank Accession Number: BC031999	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 3283	Recommended Dilutions: WB: 1:5000-1:50000
Source: Rabbit	UNIPROT ID: P14060	
Isotype: IgG	Full Name: hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 1	
Immunogen Catalog Number: AG38103	Calculated MW: 42 kDa	
	Observed MW: 42 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: human	

Background Information

HSD3B1, also known as hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 1, is an enzyme that plays a crucial role in steroid biosynthesis. It belongs to the short-chain dehydrogenase/reductase (SDR) family and is primarily involved in the conversion of C19 and C21 steroids with a 3 β -hydroxy-5-ene structure to 3-oxo-4-ene products. This enzyme is mainly localized to the placenta and non-steroidogenic tissues. It is also involved in various biological processes, including C21-steroid hormone metabolic processes, hippocampus development, and response to corticosterone. HSD3B1 has been implicated in conditions such as hypertension and hypospadias. In addition, HSD3B1 has been studied in the context of prostate cancer, where it is involved in the synthesis of androgens within the tumor, contributing to resistance to radiotherapy and other treatments. Its expression and activity can be influenced by genetic variations, which may affect clinical outcomes in patients with prostate cancer. HSD3B1 is also involved in the metabolism of oxysterols, which are oxidized forms of cholesterol or its precursors. It catalyzes the oxidation of the 3 β -hydroxy group to a 3-ketone and isomerizes the double bond from Δ 5 to Δ 4, a key reaction in bile acid biosynthesis. This enzyme's activity is essential for the conversion of initial 3 β -hydroxy stereochemistry to the 3 α -hydroxy stereochemistry in primary bile acids. HSD3B1 is also involved in the metabolism of side-chain oxysterols, which are ligands to various receptors and play roles in cholesterol biosynthesis and other biological processes. Overall, HSD3B1 is a multifunctional enzyme with significant roles in steroid hormone production, oxysterol metabolism, and various disease processes.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32642-1-AP (HSD3B1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.