

For Research Use Only

EID2 Polyclonal antibody, PBS Only

Catalog Number:32693-1-PBS



Basic Information

Catalog Number:

32693-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38433

GenBank Accession Number:

NM_153232.3

GeneID (NCBI):

163126

UNIPROT ID:

Q8N6I1

Full Name:

EP300 interacting inhibitor of
differentiation 2

Calculated MW:

25kDa 236aa

Observed MW:

28 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-EID2 (IP:32693-1-AP, 4ug; Detection:32693-1-AP 1:500) with HEK-293 cells lysate 1800 ug. This data was developed using the same antibody clone with 32693-1-PBS in a different storage buffer formulation.