

For Research Use Only

FAM171B Polyclonal antibody

Catalog Number: 32704-1-AP



Basic Information

Catalog Number: 32704-1-AP	GenBank Accession Number: BC060872	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 165215	
Source: Rabbit	UNIPROT ID: Q6P995	
Isotype: IgG	Full Name: family with sequence similarity 171, member B	
Immunogen Catalog Number: AG37311	Observed MW: 55 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

FAM171B belongs to the FAM171 family. It is an evolutionarily conserved protein with a wide range of expression across various tissues. FAM171B may play a role in stimulating immune infiltration and immune responses in diseases such as pulmonary arterial hypertension.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data