

For Research Use Only

# TRMT9B Polyclonal antibody, PBS Only

Catalog Number: 32716-1-PBS



## Basic Information

<b>Catalog Number:</b> 32716-1-PBS	<b>GenBank Accession Number:</b> BC016633	<b>Purification Method:</b> Antigen affinity Purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 57604	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9P272	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 8 open reading frame 79	
<b>Immunogen Catalog Number:</b> AG36370	<b>Observed MW:</b> 45 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
human

## Background Information

TRMT9B is one of two metazoan homologs of yeast tRNA methyltransferase 9 (Tm9), which methylates tRNA wobble uridines (PMID: 37642556).

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32716-1-AP (C8orf79 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32716-1-PBS in a different storage buffer formulation.