

For Research Use Only

BMP15 Polyclonal antibody

Catalog Number: 32721-1-AP



Basic Information

Catalog Number:

32721-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37480

GenBank Accession Number:

NM_009757

GeneID (NCBI):

12155

UNIPROT ID:

Q9Z0L4

Full Name:

bone morphogenetic protein 15

Calculated MW:

45 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 32721-1-AP (BMP15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.