

For Research Use Only

FAM3D Polyclonal antibody

Catalog Number: 32785-1-AP



Basic Information

Catalog Number: 32785-1-AP	GenBank Accession Number: NM_146050.2	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 600 ug/ml by Nanodrop;	GeneID (NCBI): 18300	Recommended Dilutions: WB: 1:1000-1:8000 IHC: 1:400-1:1600
Source: Rabbit	UNIPROT ID: P97805	
Isotype: IgG	Full Name: oncoprotein induced transcript 1	
Immunogen Catalog Number: EG2694	Calculated MW: 25 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications: WB, IHC, ELISA	Positive Controls:
Species Specificity: mouse, rat	WB : mouse colon tissue, rat colon tissue
	IHC : mouse colon tissue,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

FAM3D (family with sequence similarity 3, member D) is a protein that belongs to the FAM3 family, which is a novel cytokine-like gene family identified in 2002. This family comprises four members: FAM3A, FAM3B, FAM3C, and FAM3D. FAM3D is a gut-secreted protein that is regulated by nutritional status. It has a predicted secondary structure of a four-helix bundle, which is a common feature in many other cytokines.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32785-1-AP (FAM3D antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 32785-1-AP (FAM3D antibody) at dilution of 1:800 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 32785-1-AP (FAM3D antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).