

For Research Use Only

TAFA1/FAM19A1 Polyclonal antibody, PBS Only

Catalog Number: 32791-1-PBS



Basic Information

Catalog Number:

32791-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37253

GenBank Accession Number:

BC025746

GeneID (NCBI):

407738

UNIPROT ID:

Q7Z5A9

Full Name:

family with sequence similarity 19
(chemokine (C-C motif)-like),
member A1

Observed MW:

15 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32791-1-AP (TAF1/FAM19A1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32791-1-PBS in a different storage buffer formulation.