

For Research Use Only

RBM34 Polyclonal antibody

Catalog Number: 32834-1-AP



Basic Information

Catalog Number:

32834-1-AP

Size:

150ul , Concentration: 700 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37711

GenBank Accession Number:

BC029451

GeneID (NCBI):

23029

UNIPROT ID:

P42696

Full Name:

RNA binding motif protein 34

Calculated MW:

430 aa, 49 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HeLa cells, Jurkat cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32834-1-AP (RBM34 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.