

For Research Use Only

C19orf47 Polyclonal antibody

Catalog Number: 32841-1-AP



Basic Information

Catalog Number:

32841-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36226

GenBank Accession Number:

BC023623

GeneID (NCBI):

126526

UNIPROT ID:

Q8N9M1

Full Name:

chromosome 19 open reading frame

47

Observed MW:

45 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, Mouse skeletal muscle

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 32841-1-AP (C19orf47 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.