

For Research Use Only

C19orf47 Polyclonal antibody, PBS Only

Catalog Number: 32841-1-PBS



Basic Information

Catalog Number:

32841-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36226

GenBank Accession Number:

BC023623

GeneID (NCBI):

126526

UNIPROT ID:

Q8N9M1

Full Name:

chromosome 19 open reading frame

47

Observed MW:

45 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 32841-1-AP (C19orf47 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32841-1-PBS in a different storage buffer formulation.