

For Research Use Only

APOBEC1 Polyclonal antibody

Catalog Number: 32857-1-AP



Basic Information

Catalog Number: 32857-1-AP	GenBank Accession Number: NM_001304566	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 350 ug/ml by Nanodrop;	GeneID (NCBI): 339	Recommended Dilutions: WB: 1:500-1:3000
Source: Rabbit	UNIPROT ID: P41238	
Isotype: IgG	Full Name: apolipoprotein B mRNA editing enzyme, catalytic polypeptide 1	
Immunogen Catalog Number: AG26715	Calculated MW: 28 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, L02 cells
Species Specificity: human	

Background Information

Apolipoprotein B mRNA-editing enzyme catalytic subunit 1 (APOBEC1) is a member of the cytidine deaminase enzyme family. It forms a multiple-protein editing holoenzyme with APOBEC1 complementation factor (ACF) and APOBEC1 stimulating protein (ASP) (PMID: 30844405). This holoenzyme is involved in the editing of C-to-U nucleotide bases in apolipoprotein B and neurofibromatosis-1 mRNAs (PMID: 24916387; 11727199). APOBEC1 plays a crucial role in RNA editing, which is essential for the proper function of these proteins. Additionally, APOBEC1 has been implicated in various cellular processes, including the regulation of gene expression and the maintenance of genomic stability (PMID: 33094286).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32857-1-AP (APOBEC1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.