

For Research Use Only

HCRTR2 Polyclonal antibody, PBS Only

Catalog Number:32859-1-PBS



Basic Information

Catalog Number:

32859-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38128

GenBank Accession Number:

NM_001526

GeneID (NCBI):

3062

UNIPROT ID:

O43614

Full Name:

hypocretin (orexin) receptor 2

Calculated MW:

440aa, 51 kDa

Observed MW:

70-73 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32859-1-AP (HCRT2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32859-1-PBS in a different storage buffer formulation.