

For Research Use Only

Angiotensinogen/AGT Polyclonal antibody

Catalog Number: 32895-1-AP



Basic Information

Catalog Number:

32895-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3273

GenBank Accession Number:

NM_134432.2

GeneID (NCBI):

24179

UNIPROT ID:

P01015

Full Name:

angiotensinogen (serpin peptidase inhibitor, clade A, member 8)

Calculated MW:

52 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

rat

Positive Controls:

WB : rat liver tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32895-1-AP (Angiotensinogen/AGT antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.